



CHANGE X PLO RER

presents

ChangeClear™

A paradigm shift in technology for cell culture clarification and purification!

QUICK

Multiproduct Formats

Low Validation Costs

Quick Product
Change Over

FLEXIBLE

High Cell Titer

High Cell Viability

Any Bioreactor Scale

REMOVES

Turbidity

Endotoxins

HCP /Viruses

DNA

A cost-effective single-use technology based on Celpure® powder technology!

How to Apply Techniques to New Disposable Processes

Results were published on a poster presented by Merck Serono Biodevelopment, Martillac, France.

Introduction

Diatomaceous earth (*Baccillariophyta*)

- Silica skeleton of unicellular algae treated heat (above 1000°C), acid solutions, calcination.
- One of the components of depth filters used for cell removal and nanofilter prefiltration.
- Single Use and inexpensive technology.



Celpure®


- Highly purified diatomaceous earth with low endotoxin and low impurities level.
- GMP grade Diatomaceous earth.
- Can be directly provided in disposable bags.

	Celpure® 300 (high purity grade)	Celline Standard Super Cells (conventional grade)
Al	nd	1000
Ca	nd	1050
Mg	124	1010
Fe	56	400
Zn	nd	210
Cu	12	18
Si	nd	12
Mn	4	14
Cr	nd	4
Cd	nd	0.5

Table comparing impurities level

Results

Use of Celpure® for clarification step




- Diatomaceous earth has at least 50% more filtrability capacity than classical depth filter.
- Cost is reduced by 2 fold by using Celpure® instead of depth filter.
- Use of slurry in disposable bag (ChangeClear) facilitates Celpure® handlings.

Use of Celpure® for host cell proteins removal

- HCP removal even with cells

Protein	Filter	Without cells	With cells
Basic pI Protein	AEX	98%	97%
	AEEX	52%	97%
Acidic pI Protein	AEX	99%	67%
	AEEX	99%	ND





- HCP and DNA removal improved without cells

Discussion

Interactions with Celpure®

- Incompressible porous particles keep flow while retaining particles in the inner surface.
- Different silica forms depending on pH :
 - Ionic interactions (a), (b)
 - Hydrophobic interactions (c)
- Celpure® could replace a chromatographic resin or a charged membrane (cost saving).



Conclusions


By improving cell growth and productivity, DSP is foreseen as the future bottleneck for biotechnology companies. Whereas some industries are investing more and more in costly new technologies, Merck Serono Biodevelopment has decided to explore other inexpensive, disposable and well-known techniques used in other industries. Diatomaceous earth is one of the components of single-use depth filters used in biotechnology companies for cell clarification, nanofilter prefiltration and recently, as some studies have shown, for removal of contaminants. Highly Purified diatomaceous earth (Celpure®), with low endotoxin and low impurity levels, is now available for GMP manufacturing purposes. Cell clarification studies demonstrated that Celpure® is 2 up to 3 fold more efficient than depth filters in removing cells and can replace centrifugation step. Besides a partial impurities removal was also obtained even with cells. Moreover, the first trials on clarified harvests demonstrated the ability of such technology to remove up to 67% of impurities in the filtrate. Celpure® can be used at different purification steps and even replace classical chromatographic resins.

ChangeXplorer Production SAS, 1 Chemin du Petit Parc, 02220 BRAINE, FRANCE

Tel : +33 (0) 679 046 266 • Fax : +33 (0) 368 343 020

Email : changexplorer@orange.fr • Web : www.changexplorer.com

ChangeClear Principle Configuration



ChangeXplorer Production SAS, 1 Chemin du Petit Parc, 02220 BRAINE, FRANCE

Tel : +33 (0) 679 046 266 • Fax : +33 (0) 368 343 020

Email : changexplorer@orange.fr • Web : www.changexplorer.com

Large-scale Production System

Large-scale production is supported with a unique valve that, while disposable, prevents operator contamination.

Large-scale ChangeClear™ saves time and money by safely dispensing DE, with easy disposal, while providing a technological edge.

